


Anti-Cd355/Crtam Rabbit Monoclonal Antibody For Wb, If/Icc, Fc, Elisa - O95727

Item Number: CM0007428



Introduction

CamelBio's CD355/CRTAM rabbit monoclonal antibody validated for Western blot, IF/ICC, flow cytometry, and ELISA. Reactive with human CRTAM, a cell adhesion molecule regulating T-cell and NK cell immunity. Ideal for immunology research.

[Learn More](#)

Parameter	Value
Product Name	CD355/CRTAM Rabbit mAb
Remarks/Alias	CD355; CD355/CRTAM
Species	Human
Gene ID (Human)	56253
Gene ID	56253
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-287 of human CD355/CRTAM (NP_062550.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, FC, ELISA
Cross Reactivity	Human
SWISS	O95727
Protein Weight	22 kDa/44 kDa
Shipping	Ice bag