

Anti-Cd35/Cr1 Monoclonal Antibody (Rabbit Mab) For If, Ihc, Elisa - P17927

Item Number: CM0004104



Introduction

Rabbit monoclonal antibody targeting human CD35/CR1 (P17927). Validated for IF-P, IHC-P, mIHC, and ELISA. Useful for complement receptor studies, immune complex clearance, and viral entry research.

[Learn More](#)

Parameter	Value
Product Name	CD35/CR1 Rabbit mAb
Remarks/Alias	KN; C3BR; C4BR; CD35; CD35/CR1
Species	Human
Gene ID	1378
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 1940-2039 of human CD35/CR1 (P17927)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF-P, IHC-P, mIHC, ELISA
Cross Reactivity	Human
SWISS	P17927
Protein Weight	224kDa
Shipping	Ice bag