

Anti-Cd35/Cr1 Polyclonal Antibody For Wb, If/Icc, Elisa - P17927

Item Number: CM0018688



Introduction

Rabbit polyclonal antibody targeting human CD35/CR1 (C3b/C4b receptor, P17927).
Validated for Western blot, IF/ICC, and ELISA.
Useful for complement pathway, immune clearance, and viral entry studies.

[Learn More](#)

Parameter	Value
Product Name	CD35/CR1 Rabbit pAb
Remarks/Alias	KN; C3BR; C4BR; CD35; CD35/CR1
Species	Human
GeneID (Human)	1378
GeneID	1378
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 420-550 of human CD35/CD35/CR1 (NP_000642.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P17927
Protein Weight	224kDa
Shipping	Ice bag