

Anti-Cd36 Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q08857

Item Number: CM0029031

Introduction

High-quality CD36 rabbit polyclonal antibody, raised against a recombinant fusion protein, validated for WB, IF/ICC, IHC-P, and ELISA. Reacts with mouse and rat CD36. Ideal for studies in fatty acid metabolism, inflammation, and angiogenesis. SWISS: Q08857.

[Learn More](#)

Parameter	Value
Product Name	CD36 Rabbit pAb
Remarks/Alias	FAT; GPIV; Scarb3
Species	Mouse
GeneID (Human)	948
GeneID	12491
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-439 of mouse CD36 (NP_031669.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q08857
Protein Weight	53 kDa
Shipping	Ice bag