

Anti-Cd36/Sr-B3 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P16671

Item Number: CM0005917



Introduction

Anti-CD36/SR-B3 rabbit monoclonal antibody for research on fatty acid transport, inflammatory response, and angiogenesis. Validated for WB, IF-P, IHC-P, ELISA; cross-reacts with Human, Mouse, Rat.

[Learn More](#)

Parameter	Value
Product Name	CD36/SR-B3 Rabbit mAb
Remarks/Alias	FAT; GP4; GP3B; GPIV; CHD57; PASIV; SCARB3; BDPLT10
Species	Human
Gene ID (Human)	948
Gene ID	948
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 30-439 of human CD36/SR-B3 (NP_000063.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16671
Protein Weight	53kDa
Shipping	Ice bag