


Anti-Cd38 Monoclonal Antibody () For Fc - P28907

Item Number: CM0008039



Introduction

Rabbit monoclonal anti-human CD38 antibody conjugated with , validated for flow cytometry. Targets the multifunctional transmembrane glycoprotein CD38 (ADPRC1/cADPR1), involved in calcium signaling and immune cell adhesion. Suitable for human samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD38 mAb
Remarks/Alias	ADPRC1; cADPR1; ADPRC 1
Species	Human
GeneID (Human)	952
GeneID	952
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 43-300 of human CD38 (NP_001766.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P28907
Protein Weight	34kDa
Shipping	Ice bag