


Anti-Cd38 Monoclonal Antibody For Wb, Ihc-P, Fc, Elisa - P28907

Item Number: CM0005077



Introduction

Rabbit monoclonal antibody against human CD38, validated for WB, IHC-P, FC, and ELISA. Recognizes human and cynomolgus monkey CD38, a multifunctional enzyme and receptor key to calcium signaling and immune cell migration.

[Learn More](#)

Parameter	Value
Product Name	CD38 Rabbit mAb
Remarks/Alias	ADPRC1; cADPR1; ADPRC 1
Species	Human
Gene ID (Human)	952
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, FC, ELISA
Cross Reactivity	Human, Cynomolgus monkey
SWISS	P28907
Protein Weight	14 kDa/34 kDa
Shipping	Ice bag