

Anti-Cd38 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P28907

Item Number: CM0019898



Introduction

CD38 Rabbit Polyclonal Antibody (pAb) validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with Human, Mouse, Rat. Targets ADP-ribosyl cyclase/cyclic ADP-ribose hydrolase 1 for immunology and metabolism research.

[Learn More](#)

Parameter	Value
Product Name	CD38 Rabbit pAb
Remarks/Alias	ADPRC1; CADPR1; ADPRC 1; CD38
Species	Human
Gene ID (Human)	952
Gene ID	952
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 181-300 of human CD38 (NP_001766.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28907
Protein Weight	34kDa
Shipping	Ice bag