

Anti-Cd39/Entpd1 Monoclonal Antibody For Fc, Elisa - P49961

Item Number: CM0009321

Introduction

Recombinant rabbit monoclonal antibody targeting human CD39/ENTPD1, an ectonucleotidase that hydrolyzes ATP/ADP. Suitable for flow cytometry (FC) and ELISA, supporting research in thrombosis and immunity.

[Learn More](#)

Parameter	Value
Product Name	CD39/ENTPD1 Rabbit mAb
Remarks/Alias	CD39; SPG64; ATPDase; NTPDase-1; CD39/ENTPD1
Species	Human
Gene ID (Human)	953
Gene ID	953
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 255-442 of human CD39/ENTPD1 (NP_001767.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC, ELISA
Cross Reactivity	Human
SWISS	P49961
Protein Weight	58kDa
Shipping	Ice bag