

Anti-Cd3D Polyclonal Antibody For Wb, Elisa - P04235

Item Number: CM0023415

Introduction

Rabbit polyclonal antibody targeting mouse CD3d, validated for Western blot and ELISA. Suitable for studying T-cell receptor signaling and adaptive immunity in mouse models.

[Learn More](#)

Parameter	Value
Product Name	CD3d Rabbit pAb
Remarks/Alias	T3d; CD3d
Species	Mouse
Gene ID	12500
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 15-115 of mouse CD3d (NP_038515.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P04235
Protein Weight	19KDa
Shipping	Ice bag