

Anti-Cd3D Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P04234

Item Number: CM0018477



Introduction

Rabbit polyclonal antibody targeting human CD3D, validated for WB, IHC-P, IF/ICC, and ELISA. Detects CD3D in human and rat samples. Suitable for T-cell signaling research and diagnostic immunoassay development.

[Learn More](#)

Parameter	Specification
Product Name	CD3D Rabbit pAb
Remarks/Alias	T3D, IMD19, CD3DELTA, CD3-DELTA, CD3D
Species	Human
GeneID (Human)	915
GeneID	915
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-105 of human CD3D (NP_000723.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	P04234
Protein Weight	19kDa
Shipping	Ice bag