

## Anti-Cd3E Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P07766

Item Number: CM0013011



### Introduction

Rabbit polyclonal antibody targeting human CD3E, validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat CD3E. Suitable for T-cell receptor signaling, adaptive immunity, and immune-oncology research.

[Learn More](#)

| Parameter        | Value                                                                                                                     |
|------------------|---------------------------------------------------------------------------------------------------------------------------|
| Product Name     | CD3E Antigen Rabbit pAb                                                                                                   |
| Remarks/Alias    | T3E; TCRE; IMD18; CD3epsilon; CD3E Antigen                                                                                |
| Species          | Human                                                                                                                     |
| GeneID (human)   | 916                                                                                                                       |
| GeneID           | 916                                                                                                                       |
| Immunogen        | Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 111-207 of human CD3E (NP_000724.1). |
| Source           | Rabbit                                                                                                                    |
| Category         | Polyclonal Antibodies                                                                                                     |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                                  |
| Cross Reactivity | Human, Mouse, Rat                                                                                                         |
| SWISS            | P07766                                                                                                                    |
| Protein Weight   | 23kDa                                                                                                                     |
| Shipping         | Ice bag                                                                                                                   |