

Anti-Cd3E Rabbit Polyclonal Antibody For Wb, Ip, If/Icc, If-P, Ihc-P, Fc, Elisa - P07766

Item Number: CM0021362



Introduction

CD3E Rabbit pAb is a polyclonal antibody raised against human CD3 epsilon chain (aa 1-126). Suitable for WB, IP, IF/ICC, IF-P, IHC-P, FC, and ELISA. Cross-reacts with Human, Mouse, and Rat. Ideal for T-cell activation and adaptive immunity studies.

[Learn More](#)

Parameter	Value
Product Name	CD3E Rabbit pAb
Remarks/Alias	T3E; TCRE; IMD18; CD3epsilon; CD3E
Species	Human
GeneID (human)	916
GeneID	916
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-126 of human CD3E (NP_000724.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, IF/ICC, IF-P, IHC-P, FC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07766
Protein Weight	23kDa
Shipping	Ice bag