

## Anti-Cd3G Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P09693

Item Number: CM0010022



### Introduction

High-quality rabbit monoclonal antibody targeting human CD3G (CD3 gamma chain). Validated for Western blot, IHC-P, IF/ICC, and ELISA. Ideal for T-cell receptor signaling, adaptive immunity, and T-cell activation research.

[Learn More](#)

| Parameter        | Value                                                                                                                                  |
|------------------|----------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | CD3G Rabbit mAb                                                                                                                        |
| Remarks/Alias    | T3G; IMD17; CD3GAMMA; CD3-GAMMA; CD3G                                                                                                  |
| Species          | Human                                                                                                                                  |
| Gene ID (Human)  | 917                                                                                                                                    |
| Gene ID          | 917                                                                                                                                    |
| Immunogen        | Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 23-116 of human CD3G (NP_000064.1)) |
| Source           | Rabbit                                                                                                                                 |
| Category         | Monoclonal Antibodies                                                                                                                  |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                                               |
| Cross Reactivity | Human                                                                                                                                  |
| SWISS            | P09693                                                                                                                                 |
| Protein Weight   | 20kDa                                                                                                                                  |
| Shipping         | Ice bag                                                                                                                                |