

Anti-Cd3G Rabbit Monoclonal Antibody For Wb - P09693

Item Number: CM0006237



Introduction

CD3G Rabbit mAb is a rabbit monoclonal antibody targeting human CD3G (P09693). Suitable for WB, IHC-P, IF/ICC, and ELISA applications, with reported cross-reactivity in human and mouse. Ideal for T-cell receptor signaling and immune function research.

[Learn More](#)

Parameter	Value
Product Name	CD3G Rabbit mAb
Remarks/Alias	T3G; IMD17; CD3GAMMA; CD3-GAMMA; CD3G
Species	Human
GeneID (Human)	917
GeneID	917
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P09693
Protein Weight	20 kDa
Shipping	Ice bag