

Anti-Cd3G Rabbit Polyclonal Antibody For Wb, Elisa - P09693

Item Number: CM0019561



Introduction

High-quality rabbit polyclonal antibody against human CD3G, suitable for Western blot (WB) and ELISA. Recognizes the T-cell surface glycoprotein CD3 gamma chain, a key component of the TCR-CD3 complex involved in adaptive immunity.

[Learn More](#)

Parameter	Value
Product Name	CD3G Rabbit pAb
Remarks/Alias	T3G; IMD17; CD3GAMMA; CD3-GAMMA; CD3G
Species	Human
Gene ID (Human)	917
Gene ID	917
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-116 of human CD3G (NP_000064.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P09693
Protein Weight	20kDa
Shipping	Ice bag