

Anti-Cd3H (Cd247) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P20963

Item Number: CM0002047



Introduction

CD3H Rabbit Monoclonal Antibody validated for WB, IHC-P, ELISA in human, mouse, rat. Targets T-cell surface glycoprotein CD3 zeta chain (CD247). Suitable for T-cell signaling and adaptive immunity research.

[Learn More](#)

Parameter	Value
Product Name	CD3H Rabbit mAb
Remarks/Alias	T3Z; CD3H; CD3Q; CD3Z; TCRZ; IMD25; CD3ZETA; CD3-ZETA
Species	Human
GeneID (Human)	919
GeneID	919
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-164 of human CD3H (NP_932170.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20963
Protein Weight	19kDa
Shipping	Ice bag