

Anti-Cd3H Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P20963

Item Number: CM0024573



Introduction

Rabbit polyclonal antibody targeting human CD3H (CD3 zeta, CD247) for Western blot, immunofluorescence, immunohistochemistry, and ELISA. Cross-reacts with mouse and rat. Binds the TCR-CD3 signaling subunit.

[Learn More](#)

Parameter	Value
Product Name	CD3H Rabbit pAb
Remarks/Alias	T3Z; CD3H; CD3Q; CD3Z; TCRZ; IMD25; CD3ZETA; CD3-ZETA
Species	Human
GeneID (human)	919
GeneID	919
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 50-164 of human CD3H (NP_932170.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20963
Protein Weight	19 kDa
Shipping	Ice bag