

Anti-Cd4 Monoclonal Antibody () For Flow Cytometry - P01730

Item Number: CM0005343



Introduction

High-quality rabbit monoclonal antibody against human CD4, conjugated to fluorophore, validated for flow cytometry (FC) to detect T-helper cells and other CD4+ immune cells.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD4 mAb
Remarks/Alias	T4; IMD79; Leu-3; OKT4D; CD4mut
Species	Human
GeneID (Human)	920
GeneID	920
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P01730
Protein Weight	51kDa
Shipping	Ice bag