

Anti-Cd4 Rabbit Polyclonal Antibody For Wb, If/Icc, If-P, Elisa - P01730

Item Number: CM0013529



Introduction

Rabbit polyclonal antibody against human CD4, targeting amino acids 26-200. Validated for WB, IF/ICC, IF-P, ELISA. Cross-reacts with human, mouse, rat. Ideal for studying T-cell immune response and HIV interaction. Uniprot: P01730.

[Learn More](#)

Parameter	Value
Product Name	CD4 Rabbit pAb
Remarks/Alias	T4; IMD79; Leu-3; OKT4D; CD4mut; CD4
Species	Human
Gene ID (Human)	920
Gene ID	920
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-200 of human CD4 (NP_000607.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01730
Protein Weight	51kDa
Shipping	Ice bag