

Anti-Cd4 Rabbit Polyclonal Antibody For Wb, If/Icc, If-P, Elisa - P06332

Item Number: CM0023821

Introduction

CD4 rabbit polyclonal antibody (pAb) validated for WB, IF/ICC, IF-P, ELISA; detects human and mouse CD4. Ideal for T-helper cell, macrophage, and NK cell research, including immune response and TCR signaling studies.

[Learn More](#)

Parameter	Value
Product Name	CD4 Rabbit pAb
Remarks/Alias	L3T4; Ly-4; CD4
Species	Mouse
GeneID (Human)	920
GeneID	12504
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-394 of mouse CD4 (NP_038516.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P06332
Protein Weight	51kDa
Shipping	Ice bag