

Anti-Cd4 Rabbit Polymab® For Wb, If-P, Ihc-P, Elisa - P06332

Item Number: CM0009256

Introduction

CD4 rabbit monoclonal antibody for Western blot, immunofluorescence, immunohistochemistry, and ELISA. Detects human, mouse, cynomolgus monkey CD4. Essential for T-helper cell research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	CD4 Rabbit PolymAb®
Remarks/Alias	L3T4; Ly-4
Species	Mouse/Human
GeneID (Human)	12504/ 920
GeneID	12504
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Cynomolgus monkey
SWISS	P06332
Protein Weight	51kDa
Shipping	Ice bag