


Anti-Cd40 Polyclonal Antibody For Wb, Ihc-P, Elisa - P25942

Item Number: CM0015792



Introduction

Rabbit polyclonal antibody against human CD40 (TNFRSF5). Validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Suitable for B-cell, immune signaling, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	CD40 Rabbit pAb
Remarks/Alias	p50; Bp50; CDW40; TNFRSF5; CD40
Species	Human
GeneID (human)	958
GeneID	958
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-104 of CD40 (NP_001241.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P25942
Protein Weight	31kDa
Shipping	Ice bag