


Anti-Cd40 Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P25942

Item Number: CM0014477



Introduction

High-quality rabbit polyclonal antibody targeting human CD40 (TNFRSF5), validated for WB, IF-P, IHC-P, and ELISA applications. Recognizes human and mouse CD40, suitable for B-cell and carcinoma studies.

[Learn More](#)

Parameter	Specification
Product Name	CD40 Rabbit pAb
Remarks/Alias	p50; Bp50; CDW40; TNFRSF5; CD40
Species	Human
Gene ID (Human)	958
Immunogen	Recombinant fusion protein containing a sequence within amino acids 198-277 of human CD40 (NP_001241.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P25942
Protein Weight	31kDa
Shipping	Ice bag