


Anti-Cd40L Rabbit Polyclonal Antibody For Wb And Elisa - P29965

Item Number: CM0017096



Introduction

Rabbit polyclonal antibody targeting human CD40L (CD154), validated for WB and ELISA with cross-reactivity to mouse. A reliable tool for studying T-cell co-stimulation and B-cell activation in diagnostic and research applications.

[Learn More](#)

Parameter	Value
Product Name	CD40L Rabbit pAb
Remarks/Alias	IGM; IMD3; TRAP; gp39; CD154; CD40L; HIGM1; T-BAM; TNFSF5; hCD40L
Species	Human
Gene ID	959
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 113-261 of human CD40L (NP_000065.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P29965
Protein Weight	29kDa
Shipping	Ice bag