

Anti-Cd41/Itga2B Monoclonal Antibody For Wb, Ihc-P, Elisa - P08514

Item Number: CM0002153



Introduction

High-quality rabbit monoclonal antibody against human CD41/ITGA2B (ITGA2B/GPIIb), 113 kDa, Swiss-Prot P08514. Suitable for WB, IHC-P, ELISA. Ideal for platelet biology and hemostasis studies.

[Learn More](#)

Parameter	Value
Product Name	CD41/ITGA2B Rabbit mAb
Remarks/Alias	GT; GT1; GTA; CD41; GP2B; HPA3; CD41B; GPIIb; BDPLT2; BDPLT16; PPP1R93; CD41/ITGA2B
Species	Human
GeneID (human)	3674
GeneID	3674
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CD41/ITGA2B (P08514).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P08514
Protein Weight	113 kDa
Shipping	Ice bag