

Anti-Cd41/Itga2B Polyclonal Antibody For If/Icc, Elisa - P08514

Item Number: CM0018340



Introduction

CD41/ITGA2B rabbit polyclonal antibody validated for IF/ICC and ELISA. Recognizes human Integrin alpha-IIb, a platelet membrane receptor involved in aggregation. Ideal for thrombosis and hematology research.

[Learn More](#)

Parameter	Value
Product Name	CD41/ITGA2B Rabbit pAb
Remarks/Alias	GT; GT1; GTA; CD41; GP2B; HPA3; CD41B; GPIIb; BDPLT2; BDPLT16; PPP1R93; CD41/ITGA2B
Species	Human
Gene ID (Human)	3674
Gene ID	3674
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 160-260 of human CD41/ITGA2B (NP_000410.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P08514
Protein Weight	113kDa
Shipping	Ice bag