

Anti-Cd42A Polyclonal Antibody For Wb, Elisa - P14770

Item Number: CM0021459



Introduction

Rabbit polyclonal antibody against human CD42a (GPIX) for WB and ELISA. Targets platelet glycoprotein IX, involved in vWF-mediated adhesion and hemostasis. Recombinant immunogen, 19kDa.

[Learn More](#)

Parameter	Value
Product Name	CD42a Rabbit pAb
Remarks/Alias	GPIX; CD42a
Species	Human
GeneID (human)	2815
GeneID	2815
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 17-140 of human CD42a (NP_000165.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P14770
Protein Weight	19kDa
Shipping	Ice bag