

Anti-Cd42D (Gpv) Rabbit Polyclonal Antibody For Wb And Elisa - P40197

Item Number: CM0026416

Introduction

This anti-CD42d rabbit pAb detects human and mouse GPV for WB and ELISA. GPV is a platelet receptor mediating vWF-dependent adhesion, critical for hemostasis research and IVD development.

[Learn More](#)

Parameter	Value
Product Name	CD42d Rabbit pAb
Remarks/Alias	GPV; CD42d
Species	Human
Gene ID (Human)	2814
Gene ID	2814
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 17-477 of human CD42d (NP_004479.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P40197
Protein Weight	61kDa
Shipping	Ice bag