


Anti-Cd44 (Ko Validated) Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Fc, Elisa - P16070

Item Number: CM0002028



Introduction

KO validated anti-CD44 rabbit monoclonal antibody for human samples. Ideal for WB, IHC-P, IF/ICC, FC, and ELISA applications. Recognizes CD44 antigen, a key receptor in cell adhesion, migration, and lymphocyte homing. Recombinant fusion protein immunogen.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CD44 Rabbit mAb
Remarks/Alias	IN; LHR; MC56; MDU2; MDU3; MIC4; Pgp1; CDW44; CSPG8; H-CAM; HCELL; ECM-III; HUTCH-1; HUTCH-I; ECMR-III; Hermes-1; CD44
Species	Human
GeneID (Human)	960
GeneID	960
Immunogen	Recombinant protein; recombinant fusion protein containing a sequence corresponding to amino acids 20-178 of human CD44 (NP_000601.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, FC, ELISA
Cross Reactivity	Human
SWISS	P16070
Protein Weight	82kDa
Shipping	Ice bag