


Anti-Cd44 Monoclonal Antibody For Flow Cytometry - P15379

Item Number: CM0011675



Introduction

-conjugated rabbit monoclonal antibody against mouse CD44 (Ly-24/Pgp-1). Suitable for flow cytometry. Reacts with mouse. Target molecular weight 85 kDa.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD44 mAb
Remarks/Alias	Ly-24; Pgp-1; HERMES
Species	Mouse
Gene ID (Human)	960
Gene ID	12505
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-176 of mouse CD44 (NP_033981.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P15379
Protein Weight	85kDa
Shipping	Ice bag