

Anti-Cd44 Polyclonal Antibody (Ko Validated) For Wb, Ihc-P, If/Icc, Elisa - P16070

Item Number: CM0014875



Introduction

High-quality rabbit polyclonal CD44 antibody (KO validated) for WB, IHC-P, IF/ICC, and ELISA. Recognizes human and mouse CD44, a key cell adhesion receptor in cancer and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CD44 Rabbit pAb
Remarks/Alias	IN; LHR; MC56; MDU2; MDU3; MIC4; Pgp1; CDW44; CSPG8; H-CAM; HCELL; ECM-III; HUTCH-1; HUTCH-I; ECMR-III; Hermes-1; 44
Species	Human
GeneID (Human)	960
GeneID	960
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-178 of human CD44 (NP_000601.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P16070
Protein Weight	82kDa
Shipping	Ice bag