


Anti-Cd44 Rabbit Monoclonal Antibody For Wb, If/Icc, Fc, Elisa - P15379

Item Number: CM0007275



Introduction

High-performance anti-CD44 rabbit monoclonal antibody validated for WB, IF/ICC, FC, and ELISA. Detects mouse CD44, a key receptor for cell adhesion, migration, and T-lymphocyte homing. Ideal for immunology and cancer research.

[Learn More](#)

Parameter	Value
Product Name	CD44 Rabbit mAb
Remarks/Alias	Ly-24; Pgp-1; HERMES; CD44
Species	Mouse
GeneID (Human)	960
GeneID	12505
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-176 of mouse CD44 (NP_033981.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, FC, ELISA
Cross Reactivity	Mouse
SWISS	P15379
Protein Weight	86kDa
Shipping	Ice bag