


Anti-Cd44 Rabbit Polymab® - P15379

Item Number: CM0010480



Introduction

Rabbit monoclonal antibody targeting mouse CD44 (Swiss P15379). Validated for WB, IF/ICC, IHC-P, and ELISA. Cross-reacts with human, mouse, rat. Recognizes 40-86 kDa glycoprotein involved in cell adhesion and lymphocyte homing.

[Learn More](#)

Parameter	Value
Product Name	CD44 Rabbit PolymAb®
Remarks/Alias	Ly-24; Pgp-1; HERMES
Species	Mouse
GeneID (Human)	960
GeneID	12505
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15379
Protein Weight	40-86 kDa
Shipping	Ice bag