

Anti-Cd45 Monoclonal Antibody For Wb, Ihc-P, Mihc, Elisa - P08575

Item Number: CM0004253



Introduction

Rabbit monoclonal antibody against human CD45 (LCA, T200), suitable for Western blot, immunohistochemistry, multiplex IHC, and ELISA. Validated for leukocyte common antigen detection in research and diagnostic use.

[Learn More](#)

Parameter	Value
Product Name	CD45 Rabbit mAb
Remarks/Alias	LCA; LY5; B220; CD45; L-CA; T200; CD45R; GP180; IMD105
Species	Human
GeneID (Human)	5788
GeneID	5788
Immunogen	Recombinant protein corresponding to a sequence within amino acids 194-390 of human CD45 (NP_002829.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, mIHC, ELISA
Cross Reactivity	Human
SWISS	P08575
Protein Weight	147kDa
Shipping	Ice bag