

Anti-Cd45 Rabbit Monoclonal Antibody For If/Icc, Ihc-P, Fc, Elisa - P06800

Item Number: CM0007276

Introduction

CamelBio offers anti-CD45 rabbit monoclonal antibody, validated for mouse samples in IF/ICC, IHC-P, FC, and ELISA applications. Recombinant immunogen targets the 145 kDa CD45 protein (P06800).

[Learn More](#)

Parameter	Value
Product Name	CD45 Rabbit mAb
Remarks/Alias	loc; B220; Cd45; L-CA; Ly-5; T200; CD45R; Lyt-4; CD45
Species	Mouse
Gene ID (Human)	5788
Gene ID	19264
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-566 of mouse CD45 (NP_001104786.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, IF-F, IHC-P, FC, ELISA
Cross Reactivity	Mouse
SWISS	P06800
Protein Weight	145KDa
Shipping	Ice bag