

Anti-Cd45 Rabbit Polyclonal Antibody For Wb, Ip, If/Icc, Elisa - P08575

Item Number: CM0016545



Introduction

Rabbit polyclonal antibody against human CD45 (LCA, B220, LY5). Suitable for WB, IP, IF/ICC, IF-P, ELISA; cross-reacts with human, mouse, rat. Ideal for T-cell activation, immunology, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	CD45 Rabbit pAb
Remarks/Alias	LCA; LY5; B220; CD45; L-CA; T200; CD45R; GP180; IMD105
Species	Human
Gene ID (Human)	5788
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 611-730 of human CD45 (NP_002829.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08575
Protein Weight	147 kDa
Shipping	Ice bag