

Anti-Cd45 Rabbit Polymab® For Wb, If-P, Ihc-P, Elisa - P08575 / P06800

Item Number: CM0011088



Introduction

CD45 Rabbit PolymAb® is a recombinant rabbit monoclonal antibody targeting human/mouse CD45 (P08575/P06800). Validated for WB, IF-P, IHC-P, ELISA, it supports leukocyte biology and immunodiagnostics research.

[Learn More](#)

Parameter	Value
Product Name	CD45 Rabbit PolymAb®
Remarks/Alias	LCA; LY5; B220; CD45; L-CA; T200; CD45R; GP180; IMD105
Species	Human/Mouse
GeneID (Human)	5788
GeneID (Mouse)	19264
Immunogen	Recombinant protein of human/mouse CD45
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P08575/P06800
Protein Weight	180-240kDa
Shipping	Ice bag