

Anti-Cd45Ra Rabbit Monoclonal Antibody For Ihc-P, Elisa - P08575-8

Item Number: CM0006825

Introduction

High-quality rabbit monoclonal antibody targeting human CD45RA (isoform 6, PTPRC). Validated for IHC-P and ELISA applications; suitable for leukocyte antigen research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	CD45RA Rabbit mAb
Remarks/Alias	LCA; LY5; B220; CD45; L-CA; T200; CD45R; GP180; IMD105; CD45RA
Species	Human
Gene ID (Human)	5788
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-482 of human CD45RA (P08575-8)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	P08575-8
Protein Weight	138kDa
Shipping	Ice bag