

Anti-Cd45Rb Monoclonal Antibody For Wb - P08575-9

Item Number: CM0008294

Introduction

Rabbit monoclonal antibody against human CD45RB, validated for WB, IF/ICC, IF-P, IHC-P, FC, ELISA. Specifically detects isoform 7 of receptor-type tyrosine-protein phosphatase C (CD45RB). Ideal for flow cytometry and immunohistochemistry in diagnostic and research applications.

[Learn More](#)

Parameter	Value
Product Name	CD45RB Rabbit mAb
Remarks/Alias	LCA; LY5; B220; CD45; L-CA; T200; CD45R; GP180; IMD105; CD45RB
Species	Human
GeneID (Human)	5788
GeneID	5788
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, FC, ELISA
Cross Reactivity	Human
SWISS	P08575-9
Protein Weight	135kDa
Shipping	Ice bag