

Anti-Cd45Rb Rabbit Polyclonal Antibody For Wb, If, Elisa - P08575-9

Item Number: CM0022439

Introduction

High-quality rabbit polyclonal antibody targeting human CD45RB (isoform 7 of PTPRC), suitable for WB, IF, and ELISA, with cross-reactivity to mouse and rat. Ideal for immunodetection of leukocyte common antigen.

[Learn More](#)

Parameter	Value
Product Name	CD45RB Rabbit pAb
Remarks/Alias	LCA; LY5; B220; CD45; L-CA; T200; CD45R; GP180; IMD105; CD45RB
Species	Human
Gene ID (Human)	5788
Gene ID	5788
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-463 of human CD45RB (XP_006711537.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08575-9
Protein Weight	135kDa
Shipping	Ice bag