

Anti-Cd45Ro Rabbit Polyclonal Antibody For Wb, Elisa - P08575-4

Item Number: CM0027496

Introduction

Anti-CD45RO rabbit polyclonal antibody for WB and ELISA applications. Reacts with human CD45RO leukocyte common antigen isoform (P08575-4). Reliable tool for immunodetection research.

[Learn More](#)

Parameter	Value
Product Name	CD45RO Rabbit pAb
Remarks/Alias	LCA; LY5; B220; CD45; L-CA; T200; CD45R; GP180; IMD105
Species	Human
GeneID (Human)	5788
GeneID	5788
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-416 of human CD45RO (NP_563578.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P08575-4
Protein Weight	131kDa
Shipping	Ice bag