

Anti-Cd46 Rabbit Monoclonal Antibody For Wb, If/Icc, Fc, Elisa - P15529

Item Number: CM0011335



Introduction

Recombinant rabbit monoclonal CD46 antibody for western blot, immunofluorescence, flow cytometry and ELISA. Targets human membrane cofactor protein, a complement regulator and pathogen receptor. Also known as MCP, TLX, AHUS2, MIC10, TRA2.10.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] CD46 Rabbit mAb
Remarks/Alias	MCP; TLX; AHUS2; MIC10; TRA2.10
Species	Human
Gene ID (Human)	4179
Gene ID	4179
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-328 of human CD46 (NP_002380.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, FC, ELISA
Cross Reactivity	Human
SWISS	P15529
Protein Weight	37-44kDa
Shipping	Ice bag