

Anti-Cd46 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P15529

Item Number: CM0003042



Introduction

High-quality rabbit monoclonal antibody against human CD46 (MCP). Suitable for WB, IHC-P, ELISA. Specific for human. Useful for complement regulation, viral infection, and immune modulation research.

[Learn More](#)

Parameter	Value
Product Name	CD46 Rabbit mAb
Remarks/Alias	MCP; TLX; AHUS2; MIC10; TRA2.10; CD46
Species	Human
GeneID (Human)	4179
GeneID	4179
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CD46 (P15529).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P15529
Protein Weight	44kDa
Shipping	Ice bag