

Anti-Cd47 Rabbit Polyclonal Antibody For Wb, Elisa - Q61735

Item Number: CM0022949

Introduction

Rabbit polyclonal antibody against mouse CD47 (integrin-associated protein, IAP). Validated for Western blot and ELISA; cross-reacts with rat. Suitable for immune regulation, angiogenesis, and cardiovascular signaling studies.

[Learn More](#)

Parameter	Value
Product Name	CD47 Rabbit pAb
Remarks/Alias	IAP; Itgp; 9130415E20Rik; B430305P08Rik; CD47
Species	Mouse
GeneID (Human)	961
GeneID	16423
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-135 of mouse CD47 (NP_034711.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q61735
Protein Weight	33kDa
Shipping	Ice bag