


Anti-Cd47 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q08722

Item Number: CM0022936



Introduction

High-quality rabbit polyclonal antibody against human CD47 (IAP, OA3, MER6), validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat CD47. Suitable for cell adhesion, signal transduction, and immunoregulation studies.

[Learn More](#)

Parameter	Value
Product Name	CD47 Rabbit pAb
Remarks/Alias	IAP; OA3; MER6; CD47
Species	Human
GeneID (Human)	961
GeneID	961
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 224-323 of human CD47 (NP_001768.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q08722
Protein Weight	35kDa
Shipping	Ice bag