

Anti-Cd48 Rabbit Monoclonal Antibody For Flow Cytometry And Elisa - P18181

Item Number: CM0004799



Introduction

Rabbit monoclonal antibody targeting mouse CD48 (SLAMF2, BLAST-1), a GPI-anchored immune regulator. Validated for flow cytometry and ELISA applications. Ideal for studying T-cell and NK-cell signaling in mouse models.

[Learn More](#)

Parameter	Value
Product Name	CD48 Rabbit mAb
Remarks/Alias	BCM1; BLAST; Bcm-1; BLAST1; SLAMF2; Sgp-60; BLAST-1; MEM-102
Species	Mouse
GeneID (Human)	962
GeneID	12506
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-216 of mouse CD48 (NP_031675.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC, ELISA
Cross Reactivity	Mouse
SWISS	P18181
Protein Weight	27kDa
Shipping	Ice bag