

Anti-Cd48 Rabbit Monoclonal Antibody For Wb, Fc, Elisa - P09326

Item Number: CM0009583



Introduction

High-quality rabbit monoclonal antibody targeting human CD48 (SLAMF2). Validated for WB, FC, ELISA with cross-reactivity to human and mouse. Ideal for immune cell signaling research.

[Learn More](#)

Parameter	Value
Product Name	CD48 Rabbit mAb
Remarks/Alias	BCM1, BLAST, hCD48, mCD48, BLAST1, SLAMF2, MEM-102, CD48
Species	Human
GeneID (Human)	962
GeneID	962
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-220 of human CD48 (NP_001769.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, FC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P09326
Protein Weight	19kDa/27kDa
Shipping	Ice bag