


Anti-Cd5 Monoclonal Antibody (Conjugate) For Fc - P06127

Item Number: CM0004159



Introduction

Rabbit anti-human CD5 monoclonal antibody (conjugate) for flow cytometry validation. Recombinant immunogen (aa 25-371), detects 54kDa protein. Specific for human CD5, useful in T-cell and B-cell studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD5 mAb
Remarks/Alias	CD5; LEU1; T1; CD5 molecule
Species	Human
GeneID (Human)	921
GeneID	921
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-371 of human CD5 (NP_055022.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P06127
Protein Weight	54kDa
Shipping	Ice bag