


Anti-Cd5 Monoclonal Antibody For Flow Cytometry - P06127

Item Number: CM0012274



Introduction

Rabbit anti-human CD5 monoclonal antibody (clone) validated for flow cytometry. Recognizes human CD5, a T-cell marker involved in immune regulation. Ideal for immunophenotyping and T-cell studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD5 mAb
Remarks/Alias	T1; LEU1
Species	Human
Gene ID (Human)	921
Gene ID	921
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-371 of human CD5 (NP_055022.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P06127
Protein Weight	55kDa
Shipping	Ice bag